NorthEast Monitoring, Inc.


Digital Recorder

International User Guide

| | |
|--|----|
| Digital Recorder International User Guide – English version..... | 2 |
| Guía Internacional del Usuario del Registrador Digital - Versión en castellano (Spanish) | 6 |
| Digitales Aufzeichnungsgerät Internationales Benutzerhandbuch – Deutsche Version (German)..... | 10 |
| Digital optager - International brugervejledning - Dansk version (Danish)..... | 14 |
| Mode d'emploi international de l'enregistreur numérique - Version française (French)..... | 18 |
| Registratore digitale Guida operativa internazionale - Versione italiana (Italian) | 23 |
| Internasjonal brukerhåndbok for digitalt registreringsapparat –norsk versjon (Norwegian)..... | 27 |
| Rejestrator cyfrowy Międzynarodowy przewodnik użytkownika – wersja polska (Polish) | 31 |
| Digitaalilentimen kansainvälinen käyttöopas – suomenkielinen versio (Finnish) | 35 |
| Цифровой регистратор Руководство пользователя — русская версия (Russian) | 39 |
| Guia do Usuário Internacional para o Gravador Digital - versão em português (Portuguese)..... | 43 |
| Dijital Kayıt Cihazı Uluslararası Kullanıcı Kılavuzu – Türkçe Sürüm (Turkish) | 48 |

IP22

advancing Holter technology

NorthEast Monitoring, Inc.

141 Parker Street, Suite 111
Maynard, MA 01754 USA
www.nemon.com

phone: [+1] 978-461-3992
fax: [+1] 978-461-5991
toll free USA: 866-346-5837
email: info@nemon.com

MediMark Europe, SARL
11 rue Emile Zola
BP 2332
F-38033 Grenoble CEDEX 2
France

Digital Recorder International User Guide – English version

For use with DR200 and DR180 Series Recorders

Use this Guide in addition to the Hook-up Guide or Manual provided with your recorder

Recorder Use

The NorthEast Monitoring, Inc. digital recorder is not intended to replace real-time telemetry monitoring for patients suspected of having life-threatening arrhythmias and is not for In Vitro diagnostic use.

The NorthEast Monitoring, Inc. digital recorder is a Holter and/or Event monitor designed to facilitate the ambulatory cardiac monitoring, on order of a physician, of those patients (including infants weighing less than 10 kg.) who may benefit from such monitoring, including but not limited to those with complaints of palpitations, syncope, chest pains, shortness of breath, or those who need to be monitored to judge their cardiac functionality, such as patients who have recently received pacemakers.

The data obtained by monitoring are not analyzed at the time of recording. After the recording is complete, the data must later be copied to Holter LX Analysis and/or compatible Event system to be analyzed.

The NorthEast Monitoring, Inc. digital recorders meet the AAMI EC-11/EC-38 standard for frequency response, amplitude accuracy and timing accuracy.

Operating Environment

- +10C to +45C (+50F to +110F)
- 10% to 95% Relative Humidity
- 700 hPa to 1060 hPa pressure
- This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.
- The recorder is rated IPX0 and is not for use in wet environments.

Transport parameters

- -40C to +70C (-40 F to +158F)
- 10% to 100% Relative Humidity
- 500 hPa to 1060 hPa pressure

Batteries (not included)

The DR181 requires two and the DR200 requires one, 1.5 volt AA batteries. Do not use batteries for more than one recording. Only the following batteries are recommended for use:

- Alkaline (MN1500 or the equivalent)
- AA Eveready Lithium L91 (DR200 use only)
- AA rechargeable NiMH (nickel metal hydride) battery with a rating of at least 1500 mAh. Recharge after every use

Battery Use and Storage

- If two batteries are required, always use batteries of the same type.
- Although battery life may last longer than a recording, batteries should not be re-used for a second patient.
- Do not leave batteries in the recorder for more than two weeks when the recorder is not in use.
- Store batteries at "room" temperatures (50-90 F) and in a dry location. Do not refrigerate or freeze.
- Batteries should be disposed of following local ordinances.
- If you use rechargeable batteries:
 - The battery recharger should be kept out of the patient environment and hook-up area
 - It is best not to leave the batteries on charge for long periods of time after the charger indicates a full charge
 - Only standard chargers that are specifically rated for use with NiMH batteries should be used such as the MAHA MH-204F or Rayovac 1-Hour charger
 - Batteries that have not been used for over two weeks should be charged before use.
 - It is best to keep pairs of batteries together so they are both discharged and charged together
 - Follow warnings and recommendations that come with the charger

Leads and Electrodes

- Compatible with standard single-use silver/silver-chloride ECG electrodes
- Compatible with standard ECO Lead wires
- Patient electrodes should be left sterile in their original packaging until use - following manufacturer's instruction - and then after use, disposed of following local ordinances and the manufacturer's instruction
- Use only NorthEast Monitoring, Inc. shielded patient cables for connecting lead wires to the recorder
- Patient lead wires and cables should be visually inspected between uses for worn or cracked areas, and
- Worn lead wires and cables should be replaced before next use and disposed of following local ordinances and the manufacturer's instructions
- Frequently used cables should be replaced at least every 6 months
- Do not remove the cable from the recorder after each use
- Do not wrap the lead wires tightly around the recorder after each use. This can damage them.

SD Cards

- The patient's data is stored in non-volatile memory internal to the recorder.
- For **Holter** recordings, patient data is stored in the recorder on a removable SD Card.
- For **Event** recordings, data is stored in non-volatile memory internal to the recorder and SD cards are not required.
- There are some SD card types that may draw excessive power, and will drain the battery prematurely. Therefore, if you purchase cards from a supplier other than NorthEast Monitoring, Inc., it is recommended that you first test the SD card for a greater amount of time than the expected use.
- If you have a new SD card that has not been formatted with a flash.dat file, you will need to use your NorthEast Monitoring, Inc. Holter LX Analysis software to initialize the card for the first time. If your software does not refer to your recorder, follow instructions for the SD360.

Maintenance and Care of Recorder, Leads and Cables

- Always remove the battery before cleaning the recorder.
- Do not leave batteries in the recorder for extended periods (more than two weeks) when the recorder is not in use.
- Clean the outside of the recorder with a damp soft cloth between uses; use water and a nonabrasive liquid soap, as required. DO NOT use any abrasive cleaners, such as acetone, on the outside of the recorder.
- Patient lead wires and cables should be cleaned and disinfected as required between uses, with 70 percent isopropyl alcohol and a soft cloth to disinfect, and a non-abrasive liquid soap and soft cloth to clean. Sterilization is not required.
- Disinfect as needed, following instructions from your infection control department. Sani-Cloth germicidal surface wipes are recommended. Sterilization is not needed.
- Do not pull on or stretch the patient cables or lead wires when you clean them. This can cause premature failure of the cable. Instead, lay the cable and attached wires on a clean, flat surface, hold them down with one hand, and holding a cloth in the other hand, rub all surfaces of the cable.
- Do not submerge the recorder or cables in water.
- Do not store the recorder at temperatures below 50 degrees F or above 110 degrees F (between 10 and 45 degrees C).
- At the end of their useful lives, all NorthEast Monitoring, Inc. products should be disposed of following local ordinances.

LCD screen and Error Messages

- **Screen shows the time:** If you see the time on the screen, the recorder is recording.
- **Battery %** -Shows the battery life remaining.
- **Blank Screen:** If the LCD screen is completely blank, this means that the recorder is not recording. It is possible that the battery no longer has enough power or batteries removed.
- **Battery LOW:** Battery is running low. When this message first appears, you will have about 5 days of recording time left on your recorder.
- **Battery FAILURE:** Recording has stopped, as the battery does not have enough power.
- **LEAD LOOSE:** This error will occur when there is a problem with the patient hook up. The problem may be with an electrode, a lead or the cable that connects the leads to the DR200. The LEAD LOOSE message will remain on the screen for about 10 seconds after the problem has been corrected. This error message can be turned off in General Settings.
- **SD Card is full:** The card cannot write any more data. Recording has stopped.
- **No card or flash.dat in recorder:** For Holter recording, an SD card with a flash.dat file is required. View card using Explorer to see files, and/ or reinitialize card using LX analysis.
- **SD Card is write locked:** Write Lock tab is set on the SD Card. Unlock Write Lock tab and try again.
- **SD Setup Failure:** Failure during write of patient ID to SD Card. You will need to re-initialize your card using the LX software.
- **SD Card Incorrectly erased:** There may be disallowed files on the SD Card. Remove SD Card from recorder and use card reader and Explorer to identify and delete these files. The only file allowed is flash.dat.
- **Unable to write to SD Card:** An error was found while attempting to write to the SD Card. This message occurs when the card is full. Sometimes this message will appear when a card is defective.

- **Write Timeout error:** This usually means a defective card.

To Remove Belt Clip (DR200 models)

If you need to remove the belt clip, you will need a long flat tool like a screw driver. In order to remove the clip, one has to slightly pry up the end of the clip near the battery cover while pulling the clip out.

DR180 Oximetry

For Oximetry recording, channels 1 and 2 are ECG, and channel 3 will be replaced with the SPO₂. If the SPO₂ cable is not attached properly, you will get the error message "No cable". If the cable is connected but without a sensor (finger-wrap, clip, etc.), you will get the error message, "No sensor".

Holter Memory Requirements

The following table lists the SD card memory requirements for recording in different modes and criteria for a 24 hour period:

| Recording mode | 12-Lead Settings: 12-lead strips saved and sampling rate | Memory required Normal Resolution 24 hours | Memory required High Resolution 24 hours |
|--------------------------|---|---|---|
| 3 channel or oximetry | n/a | 26 MB | 52 MB |
| 12-lead | 3 sec every 180 at 720 s/s | 35 MB | 61 MB |
| 12-lead | 3 sec every 18 at 720s/s | 119 MB | 145 MB |
| 12-lead | 10 sec every 60 at 720s/s | 119 MB | 145 MB |
| 12-lead | continuous at 720s/s | 585 MB | 611 MB |
| 12-lead | 3 sec every 180 at 1440s/s | 45 MB | 71 MB |
| 12-lead | 3 sec every 18 at 1440s/s | 212 MB | 238 MB |
| 12-lead | 10 sec every 60 at 1440s/s | 212 MB | 238 MB |
| 12-lead | continuous at 1440s/s | 1145 MB | 1171MB |
| SAECG | n/a | 445 MB | 471 MB |

Repairing the Recorder

Should the NorthEast Monitoring, Inc. recorder fail to work properly during its useful life, contact your local representative or distributor; or call NorthEast Monitoring, Inc. at 866-346-5837 (USA only) or [+1] 978-461-3992, option 1, to arrange for repair and/or replacement.

- There are no user-serviceable parts.
- Removal of the label to open the recorder voids the warranty.
- Prior to returning a recorder, you must obtain a Return Merchandise Authorization (RMA) number. If the RMA number is not visible on the outside of the packing carton, NorthEast Monitoring, Inc. will refuse delivery.

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- Tüm NorthEast Monitoring, Inc. ürünleri, faydalı ömürlerinin sonunda yerel yönetmeliklere göre bertaraf edilmelidir.