

NorthEast Monitoring, Inc.

advancing Holter technology™

HE/LX[®] (helix) Analysis Software - Features Matrix

Software Feature / Arrhythmia	Software Level				
	Pro	Enhanced+	Enhanced	Basic	Edit
International Languages	✓	✓	✓	✓	✓
1 Strip Per Hour	✓	✓	✓	✓	
Bi-Directional Scanning	✓	✓	✓	✓	
Diary Entry	✓	✓	✓	✓	✓
Event Marker	✓	✓	✓	✓	
Label Strips	✓	✓	✓	✓ ^{limited}	
Maximum and Minimum HR	✓	✓	✓	✓	
Disable SVPBs	✓	✓	✓	✓	
Page Mode	✓	✓	✓	✓	
Six Minute Walk Assessment	✓	✓	✓	✓	
Sleep Apnea	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	
SpO ₂	✓ ²	✓ ²	✓ ²	✓ ²	
Calibration Strips	✓	✓	✓		
Configurations	✓	✓	✓		
Critical Events	✓	✓	✓		
Failure To Capture	✓	✓	✓		
Failure To Sense	✓	✓	✓		
Inhibition	✓	✓	✓		
Insert Beat	✓	✓	✓		
Interval Editor	✓	✓	✓		
Pause	✓	✓	✓		
QTc – 4 choices for calculating	✓ ¹	✓ ¹	✓ ¹		
Saved Strips	✓	✓	✓		
ST Baseline	✓	✓	✓		
ST Depressions and Elevations	✓	✓	✓		
ST Deviation (Maximum)	✓	✓	✓		
SVPB, SVPB Pair, SVT	✓	✓	✓		
Time Domain HRV	✓	✓	✓		
True 12-Lead	✓ ¹	✓ ¹	✓ ¹		
VPB, VPB Pair, VTAC	✓	✓	✓		
AFib On/Off	✓	✓			
Receive ECG data/Transmit report files	✓	✓			
Review Bins	✓	✓			
Templates	✓	✓			
Beats	✓	✓			
EDIT mode available	✓				✓
Afib Events	✓				
Bigeminy	✓				
Bradycardia and Tachycardia	✓				
Early VPB	✓				
Fastest and Longest SVT	✓				
Fastest and Longest VTAC	✓				
Hospital Information System File	✓				
Irregular and Longest R-R	✓				
Quadrigeminy	✓				
R on T VPB	✓				
Spectral Analysis	✓				
Spectral Plots & Tables	✓				
ST End	✓				
ST Maximum HR	✓				
ST Onset	✓				
Trigeminy	✓				

- 1 - DR180 Series recorder w/ 12-Lead capability required
- 2 - DR180 Series and purchase of SpO₂ cable required
- 3 - LX Sleep option required



REPORT Features Matrix

Report Feature	Software Level				
	PRO	ENHANCED+	ENHANCED	BASIC	EDIT
Back-up and restore (CD-RW required)	✓	✓	✓	✓	
Comments page	✓	✓	✓	✓	✓
Full disclosure - full or partial	✓	✓	✓	✓	✓
Full disclosure - single or two channel	✓	✓	✓	✓	✓
Full disclosure strips	✓	✓	✓	✓	✓
Full-sized strips	✓	✓	✓	✓	✓
International languages	✓	✓	✓	✓	✓
List of diary events	✓	✓	✓	✓	✓
Oxy and HR trend - 24 hours (w/Apnea ³)	✓ ²	✓ ²	✓ ²	✓ ²	✓ ²
Oxy and HR trend - 2 hours	✓ ²	✓ ²	✓ ²	✓ ²	✓ ²
Oxy and heart rate summary	✓ ²	✓ ²	✓ ²	✓ ²	✓ ²
Oxy values with full disclosure	✓ ²	✓ ²	✓ ²	✓ ²	✓ ²
Oxy trend with full disclosure	✓ ²	✓ ²	✓ ²	✓ ²	✓ ²
Oxy respiration (Apnea patients only)	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}
Oxy airflow (Apnea patients only)	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}	✓ ^{2&3}
Patient information summary: 7 options	✓	✓	✓	4 types	✓
Review and edit front page	✓	✓	✓	✓	✓
Six-Minute-Walk Assessment	✓	✓	✓	✓	✓
Status window	✓	✓	✓	✓	✓
12-lead strips - 25 or 50 mm/sec	✓ ¹	✓ ¹	✓ ¹		✓ ¹
General profile and trends	✓	✓	✓		✓
HRV time summary	✓	✓	✓		✓
Paced data information / heart rate trend	✓	✓	✓		✓
Saved strips list	✓	✓	✓		✓
ST episodes	✓	✓	✓		✓
Time domain table	✓	✓	✓		✓
Normal bins	✓	✓			✓
Paced bins	✓	✓			✓
Ventricular bins	✓	✓			✓
12-lead table	✓ ¹				✓ ¹
12-lead trends	✓ ¹				✓ ¹
3-D plot "waterfall graph" (spectral)	✓				✓
Bigeminy Table	✓				✓
Bradycardia/Tachycardia Table	✓				✓
Circadian plots (spectral)	✓				✓
Critical events	✓				✓
Frequency domain table (spectral)	✓				✓
Hourly rhythm page	✓				✓
Normalized frequent domain (spectral)	✓				✓
Paced interval histograms	✓				✓
Paced summary information	✓				✓
Expanded HR and ST trends	✓				✓
Settings / Scanning criteria	✓				✓
Spectral HRV (frequency domain)	✓				✓
ST analysis triaxial graph	✓ ¹				✓ ¹
Supraventricular summary	✓				✓
Ventricular summary	✓				✓


- 1 - DR180 Series recorder w/ 12-Lead required
- 2 - DR180 Series w/ purchase of SpO₂ cable required
- 3 - LX Sleep option required

Industry-leading 3-year warranty

NorthEast Monitoring, Inc.

510(k)K930564

NB: In the U.S.A., Federal Law restricts devices to sale by or on the order of a physician. Not for in vitro diagnostic use. Not intended to replace real-time telemetry monitoring for patients suspected of having life-threatening arrhythmias.



advancing Holter technology

141 Parker Street, Suite 111
Maynard, MA 01754 USA
www.nemon.com

phone: [+1] 978-461-3992
fax: [+1] 978-461-5991
toll free USA: 866-346-5837
email: info@nemon.com



0413

October 2016

NEMM045-REV-C