





Compare Recorders

 nemon.com
 (+1) 978-461-3992
 Toll Free USA 866-346-5837
 sales@nemon.com

	DR400 patch style Holter recorder with event capabilities	DR300 Holter and Event Recorder with Bluetooth Wireless	DR200/HE Holter and Event Recorder	DR181 3-Channel DR181 12-Lead with OxyHolter Capability
	<p>Weighing in at just 34 grams, the showerproof 1- or 3-lead patch style recorder solves your remote cardiac monitoring needs</p> 	<p>A revolutionary recorder combining the best of Holter, Event and Mobile Cardiac Telemetry (MCT)</p> 	<p>The industry's smallest, single unit Holter and Event recorder</p> 	<p>High-quality research and clinical evaluation tool. Designed for ease-of-use, accuracy of data and durability of the unit</p> 
Current version	5.02	4.46	1.11 and 4.47	3.05
Holter Recorder	X	X	X	X
Holter Capacity	14 days	14 days	14 days	3 days
Event Recorder	X	X	X	
Event, auto-detect	X	X	X	
Event capacity	30 days	30 days	30 days	
Event Bluetooth	X	X		
Event w/ Holter back up	14 days	14 days	14 days	
Pacemaker detection	X	X	X	X
1-lead hook-up	patch electrode	X	X	
3-lead hook-up	patch electrode	X	X	
5-lead hook-up		X	X	X
7-lead hook-up		X	X	X
12-lead hook-up				optional
12-lead capacity				3 days
Oximetry				X
User interface	PCPatch utility on PC	LCD screen	LCD screen	LCD screen
Data Transfer	USB cable	SD Card	SD Card	SD Card
Industry-leading warranty	3-year	3-year	3-year	3-year
Batteries	Rechargeable internal	1- AA	1-AA	2-AA
Sampling rate	180 samples/second	180 samples/second	180 samples/second	1440 samples/second

NorthEast Monitoring is a medical electronics company that develops high quality, rugged products to address diagnostic, clinical and research needs in cardiology, respiratory, and pulmonary applications.

Contact NorthEast Monitoring for more information or to place an order.