



Compare Recorders

-  nemon.com
- (+1) 978-461-3992
- Toll Free USA 866-346-5837
- sales@nemon.com

DR400 patch style Holter recorder with event capabilities

Weighing in at just 34 grams, the showerproof 1- or 3-lead patch style recorder solves your remote cardiac monitoring needs



DR300 Holter and Event Recorder with Bluetooth Wireless

A revolutionary recorder combining the best of Holter, Event and Mobile Cardiac Telemetry (MCT)



DR200/HE Holter and Event Recorder

The industry's smallest, single unit Holter and Event recorder



Current version	5.04	4.46	1.11 and 4.47
Holter Recorder	X	X	X
Holter Capacity	14 days	14 days	14 days
Event Recorder	X	X	X
Event, auto-detect	X	X	X
Event capacity	30 days	30 days	30 days
Event Bluetooth	X	X	
Event w/ Holter back up	14 days	14 days	14 days
Pacemaker detection	X	X	X
1-lead hook-up	patch electrode	X	X
3-lead hook-up	patch electrode	X	X
5-lead hook-up		X	X
7-lead hook-up		X	X
12-lead hook-up			
12-lead capacity			
Oximetry			
User interface	PCPatch utility on PC	LCD screen	LCD screen
Data Transfer	USB cable	SD Card	SD Card
Industry-leading warranty	3-year	3-year	3-year
Batteries	Rechargeable internal	1- AA	1-AA
Sampling rate	180 samples/second	180 samples/second	180 samples/second

NorthEast Monitoring is a medical electronics company that develops high quality, rugged products to address diagnostic, clinical and research needs in cardiology, respiratory, and pulmonary applications.

Contact NorthEast Monitoring for more information or to place an order.