



DR400 patch style recorder

Revision Chart

	Rev A	Current version
		
Charging and uploading	USB into back of recorder	Docking station
Current version	n/a	5.25
Firmware	5.03 – 5.07	5.21 and newer
Current PC Patch version	1.08	1.10
Button on front says	EVENT	PRESS
When plugged into PC, LED light ...	Indicates PCPatch responses	Indicates charging status (below)
When plugged into power, LED light...	Flashes orange when charging. Flashes green when fully charged.	
Event recorder?	Yes	MCT/Event available



NorthEast Monitoring, Inc.

141 Parker Street, Suite 200
 Maynard, MA 01754 USA
www.nemon.com

phone: [+1] 978-461-3992
 fax: [+1] 978-461-5991
 toll free USA: 866-346-5837
 email: info@nemon.com

advancing Holter technology